



# Using microtomography in geosciences

Ex: Studying fluid and mass transfers in porous geomaterials with application to CO<sub>2</sub> underground storage, ...

Philippe Gouze





## Using microtomography in geosciences

Ex: Studying fluid and mass transfers in porous geomaterials with application to CO<sub>2</sub> underground storage, ...

Philippe Gouze (Transport in Porous Media research unit)

- CNRS (Governmental Research Institute): ≅ 20000 people; all disciplines
- Géosciences Montpellier :

Research covers a large range of topics in geology, geophysics and hydrogeology  $\cong$  150 people (1/3 academics, 1/3 CNRS), publications  $\cong$  130/yr, PhD defense  $\cong$  15/yr.

Based at the *University of Sciences - Montpellier* 

 $\cong$  4500 people (3/5 academics),  $\cong$  16000 students, 7 departments ( $\cong$  90 research units)

#### To make a long story short:

CO<sub>2</sub> storage (i.e. massive injection of CO<sub>2</sub> in underground reservoirs or aquifers) produces highly aggressive fluids and strong mechanical stresses.

The standard *reservoir* models developed, for instance, for the oil industry cannot handle these conditions.

Many coupled mechanisms taking place in these conditions are still poorly known.

The acceptability (both in term of industrial feasibility and risk assessment) of this technology requires sound predictive models.


#### The main issues are:

- Identifying all the mass transfer processes
   (ex: fluid-rock reactions)
- Measuring parameters for feeding the numerical models (ex: permeability, dispersivity)
- Determining the functional relationships required for coupling flow, transport and reaction in the numerical models (ex: porosity-permeability laws)

#### In this context X-ray microtomography is used to study

Media composition (rock-forming minerals and void distribution)
Media structure *versus* flow properties
Media structure *versus* solute transport & reactions
Structure of fluid-fluid interfaces (multiphase flow)
Media structure *versus* mechanical properties


media = porous media or fractured media = reservoir rocks, caprocks, well cement, ...




Reproduce in situ conditions

# fractures

#### Rough fracture dissolution: geometrical properties versus aperture increase



#### Aperture spatial correlation




#### Many other parameters ...

$$\tau = \frac{1}{l} \sum_{y=1}^{y=l} \frac{L'}{L} = \frac{1}{l} \sum_{y=1}^{y=l} \frac{\sum_{z=1}^{z=L} \sqrt{(h_{z+1} - h_z)^2 + \Delta z^2}}{L}$$


$$\Re = \sqrt{\frac{1}{Ll} \sum_{y=1}^{y=l} \sum_{z=1}^{z=L} (h_{z+1} - h_z)^2}$$

#### Average aperture // and $\perp$ to flow




(Noiriel and Gouze, JH, in review)


#### Rough fracture dissolution: flow field properties versus aperture increase




# Planar fracture (sawn) : localization effects and dissolution rate versus Q




#### Rough fracture dissolution with coating growth (marl)






(Noiriel et al., WRR, 2007)

#### Rough fracture dissolution with coating growth (marl)




(Noiriel et al., WRR, 2007)


#### Rough fracture dissolution: effect of differential kinetics




# Rough fracture dissolution: effect of differential kinetics





#### **Alteration of fractured class G well cements**





#### Alteration of fractured class G well cements














#### Alteration of fractured class G well cements




# Porous media

#### Main issue for model = poro-perm relationship






# Characterization of the dissolution/precipitation processes in the reservoir





Cumulative volume of rock dissolved at the end of the experiment



#### and measuring effective dissolution rate and kinetics parameters


## Oxfordian crinoïdal limestone (OCL)






#### and measuring effective dissolution rate and kinetics parameters


# Oxfordian crinoïdal limestone (OCL)




#### and measuring effective reactive surface area



# **Compute permeability tensor**





	$K_{ij}$	
1.216	4.77 10 <sup>-3</sup>	1.90 10 <sup>-2</sup>
4.78 10 <sup>-3</sup>	1.227	-0.173
1.90 10 <sup>-2</sup>	-0.173	1.046

#### **Compute permeability tensor**

(PC, Bernard, ICMCB, Bordeaux)


1) Meshing and then compute the Stokes equation ... still a challenge for "large" size structures: new algorithm, // computing, ...

The choice of the resolution (i.e. the sample size) is very important




 $K_{ij}$ 1.216 4.77 10<sup>-3</sup> 1.90 10<sup>-2</sup>
4.78 10<sup>-3</sup> 1.227 -0.173
1.90 10<sup>-2</sup> -0.173 1.046


# Use properties for generating statistically equivalent porous media



## **Construction of pore network models**


Objective: simplifying computations for flow and transport by using much simple geometries.







#### **Construction of pore network models**

Objective: simplifying computations for flow and transport by using much simple geometries.




# Multiphase flow ....







# Objective: measure residual saturation & relative permeability




#### Multiphase flow ....



#### Dissolution and particles transport in carbonates



pores < resolution ( $\mu$ –porosity)



Porosity always increases (i.e. solid mass is removed) but

Permeability may increase or decrease depending on the pH of the fluid.

 $\Rightarrow$  Calcite  $\mu$ -grains detached and moved instead of being dissolved


(Garing et al., subm.)


# Applications to subsurface processes (dispersion of pollutants)

#### A main issue:


Explain the always-observed non-Fickian dispersion behavior (while all models assume Fickian dispersion)

#### 1) Field scale experiment






Immobile domain



#### 3) Compare G(t) computed from the XRMT and G(t) obtained from the field measurements



#### **Main issues**

Media composition (mineral and void distribution)

Media structure versus flow properties

Media structure versus solute transport & reactions

Structure of fluid-fluid interfaces (multiphase flow)

Media structure versus mechanical properties

## Requirements (not all achievable so far, or at least not simultaneously ...)

3D!

Using samples of characteristic size around 1 cm.

Large range of scales (resolution from sub-100nm to tens of microns)

Fast imaging for allowing dynamic experiments

Resolving low contrast interfaces

**Resolving mineral composition** 


If you think about the unique capabilities of SRX

what experiment comes immediately to your mind that you would like to do? Static multi scale observations (zooming-in)

Dynamic experiments (monophasic and multiphasic fluid flow in porous media)

How should the sample environment look like?
What would be the demands on the environment from your samples?
For dynamic experiment, we need to build specific confinement cells (aluminum, carbon, ...?)

Which elements with absorption edges in the energy range covered by SRX (4.65keV to 22keV) would be of most interest for you?



If you think about the unique capabilities of SRX

what experiment comes immediately to your mind that you would like to do? Static multi scale observations (zooming-in)

Dynamic experiments (monophasic and multiphasic fluid flow in porous media)

How should the sample environment look like?
What would be the demands on the environment from your samples?
For dynamic experiment, we need to build a confinement cell (aluminum, carbon, ...?)

Which elements with absorption edges in the energy range covered by SRX (4.65keV to 22keV) would be of most interest for you?

What would be more important for you for an experiment at SRX, very high spatial resolution (sub-100nm) or a large sample area (mm with sub-micron resolution)? The ideal for us is from sub-100nm to 5 microns (using "zooming in)

